



@gzlongest  
 export@longest.cn

@gzlongest  
 @gzlongest

LGT-2320BE Specification	
Output Channel	Ⓐ 12 Channels consist by two groups; Ⓑ Each group can be set and adjusted independently.
Treatment Mode	Nine Modes Standard、Toning、Fast Slim、Lymph Drainage、Relaxation、 Abdominals、Thighs&Hips、Facial、Bust
Output Intensity	Adjustable, 0-120mA
Treatment Time	1-99min, stepping 1min
Speed (Training frequency)	Slow, Medium and Fast
Safety Protection	Over-current Protection Over-voltage Protection Short Circuit Protection



**Guangzhou Longest Science & Technology Co., Ltd.**  
**Address:** 5&6F, Building B4, No.11, Kaiyuan Avenue, Science City,  
Guangzhou Hi-tech Industrial Development Zone, 510530 Guangzhou, China  
**Tel:** 86-20-66353999   **Fax:** 86-20-66353920  
**E-mail:** export@longest.cn   **Website:** www.longest.cn

# Electronic Muscle Stimulation for Beauty




- Lose weight
- Tone your body
- Get rid of cellulite



• LGT-2320BE



**Electronic Muscle Stimulation for Beauty**  
(Model Number:LGT-2320BE) applies low-frequency currents (modulation or non-modulation) to stimulate the motor nerves or muscles to make them contract. After the muscle contractions, the muscle fiber thickens, the weight of the muscle volume increases, the capillaries in the muscle become rich, which strengthen the entire muscle to achieve the effect of lose weight and shaping.



Features

- Multiple Treatment Areas
- No any weights during your workout
- Intensity automatic locking for your safety
- Non-Invasive , No Downtime
- 45-minutes session can burn up to 1000 calories

Application

- |                   |   |
|-------------------|---|
| Back of the Chest | <b>Relaxation Mode</b><br>Build up Muscle, Strengthening and Relaxing             |
| Upper Arms        | <b>Standard Mode</b><br>Build up Muscle, Strengthening and Toning                 |
| Abdomen           | <b>Abdominals Mode</b><br>Cellulite Reduction, Strengthening and Toning           |
| Buttocks          | <b>Thigh&amp;Hip Mode</b><br>Fat Breakdown and Toning                             |
| Thigh             | <b>Thigh&amp;Hip Mode</b><br>Break down fat, Strengthening and Lymphatic Drainage |
| Lower Leg         | <b>Toning Mode</b><br>Toning and Cellulite Reduction                              |

Toxin Release

Build up Muscles

Body Reshaping

Fat Loss

Alleviate Joint Pressure

**Benefits From EMS Training**

**Before**

**Outcomes from EMS Training**

**After**

Straps with 3 Sizes\*12 (4 pcs per size)

Electrode Lead Hose (Φ95mm) \*12

**Accessories**

Silicone Rubber Electrode\*12

Safety Switch\*1